# 1. Warnings/Precautions

## WARNING

- Make sure that the personnel are made familiar with the contents of this manual. The Ambu TwinPump should only be operated by persons who have received adequate training in practising pharyngeal and tracheal suction techniques.

## CAUTION

- U.S. federal law restricts this device to sale by or on the order of a physician (USA and Canada only).
- Following each cleaning and reassembly the pump should be tested for correct function (see section 9).

# 2. Introduction

The **Ambu® TwinPump** and the **Ambu® TwinPump 1000** are high performance suction pumps designed for pharyngeal and tracheal suction in emergency situations. The pumps can be operated by foot or by hand.

The Ambu TwinPump is a double acting piston pump with high suction performance due to a combination of large airflow and high vacuum. Thanks to its special twin-cylinder design the pump is easily operated in a kneeling position, the foot resting continuously on the pump pedal without being lifted in contrast to other foot operated pumps in which the return movement is spring activated.

The compact dimensions of the Ambu TwinPump makes it the ideal pump for emergency medical cases and other locations where space is limited.

# 3. Specifications

The Ambu TwinPump and Ambu TwinPump 1000 are in conformity with the following standard: EN ISO 10079-2:1996 E

The Ambu TwinPump and Ambu TwinPump 1000 are in conformity with Council Directive 93/42/EEC concerning Medical Devices 

**Dimensions (without suction tube):** Length: 206 mm, Width: 96 mm, Height: 104 mm/146 mm (TwinPump 1000)

**Weight (incl. suction tube and tip):** Approx. 1.1 kg/1.2 kg (TwinPump 1000)

**Max. vacuum:** Approx. -80 kPa (-800 cmH<sub>2</sub>O/-600 mm Hg)

**Free airflow at two pumping strokes per second:** Approx. 30-40 litres/min.

**Free airflow at max. pumping strokes:** Approx. 70 litres/min.

**Usable capacity of collection container:** Approx. 600 ml/1000 ml (TwinPump 1000)

**Internal diameter of suction tube:** 10 mm

**Operating temperature:** -20 °C to 50 °C

**Storage temperature:** -40 °C to 70 °C

## Materials

### Parts

Transparent plastic parts  
Bottom cover  
Manifold pipe  
Gaskets, O-rings and valve diaphragms  
Piston rings  
Foot pedal  
Other metal parts  
Suction tip  
Suction tube

### Materials

Polycarbonate  
Thermoplastic rubber  
Polypropylene  
Silicone rubber  
Teflon  
Aluminium  
Nickel plated brass and stainless steel  
Acetal  
Silicone rubber

## 4. Principles of operation - see ①

The pump is operated by seesaw movement of the pedal.

Depressing the left-hand side of the pedal causes the valve diaphragm (1.1) to open, thus expelling the air from the cylinder (1.2). The valve diaphragm (1.3) is closed. At the same time, the right-hand side of the pedal rises, causing the valve diaphragm (1.4) to open and the air in the collection container (1.5) to be drawn into the cylinder (1.6). This creates a vacuum in the container (1.5). The valve diaphragm (1.7) is closed. Air and fluid therefore flow from the suction tip through the suction tube and the suction tube connector (1.8) into the container (1.5).

A corresponding process takes place when the right-hand side of the pedal is depressed.

As suction is generated every time one side of the pedal is depressed, the resulting flow is substantially greater - approx. twice the normal value - than for manual pumps, which only develop a negative pressure during one half of the working cycle.

## 5. Operating instructions

The pump can be operated by foot or by hand. See ②.

Work the pump using continuous seesaw movements on the pedal. Very little force is required to maintain a high airflow through the suction tube. When the negative pressure increases, this is immediately noticeable as resistance to pedal movement.

The collection container has a fluid capacity of approx. 600 ml. for TwinPump and 1000 ml. for TwinPump 1000. If the container becomes full at a critical moment and there is no time to empty it, do not interrupt the suction. The pump will continue to operate as the excess fluid will be discharged through the two valve diaphragms on top of the cylinders.

The suction tube is fitted with a combination suction tip. The small bore tip can be used with or without a suction catheter. See ③. To accommodate large quantities of fluid and solids, e.g. patient vomit, the small bore tip is easily removed to enable instant use of the large bore tip (10 mm aperture). This will substantially increase the rate of suction and permit the passage of solids without obstructing the aperture. See ④. During the process the small bore tip should be hooked securely to the suction tube to prevent loss.

## 6. Disassembly

For hygienic reasons the use of gloves is recommended during disassembly and cleaning.

Pull away and remove the bottom cover. See ⑤. Rinse out the contents of the collection container. To clean thoroughly, the pump may now be disassembled as shown in ⑥ - ⑧. Furthermore the piston ring and the O-ring mounted on each piston should be detached. The valve diaphragms should not be removed.

## 7. Cleaning - disinfection - sterilisation

After disassembling the pump as described in section 6, the parts can be cleaned, disinfected and sterilised as follows:

### 7.1. Manual washing

Wash the parts in hot water and detergent. Now rinse carefully and make sure that all traces of detergent have been removed. Allow the parts to dry before reassembly.

### 7.2. Machine washing

The washing machine must be designed for handling medical instruments and anaesthesia equipment. Use the washing programme and detergent normally employed for anaesthesia equipment. Allow the parts to cool and dry completely before reassembly.

**Note:** The aluminium pedal should not be cleaned by a washing machine, since the detergent will oxidize the surface.

### 7.3. Autoclaving

The parts are autoclavable at a temperature of 121 °C when autoclaved separately.

### 7.4. Boiling

Place the parts in boiling water and boil for 10 minutes. Allow cooling and drying before reassembly.

### 7.5. Chemical disinfection

Disinfectants must be suitable for the materials used in pump manufacture (see section 3). The supplier's instructions concerning mixing proportion and disinfecting time should be closely adhered to. It is important to rinse all parts carefully in clean water so as to remove all traces of disinfectant that might make parts sticky or reduce their service life. Allow the parts to dry before reassembly.

### 7.6. Gas sterilisation

Gas sterilisation should not be applied. Polycarbonate surfaces may become sticky and lock the valve diaphragms.

## 8. Assembly

Check that all pump components are undamaged.

Mount the O-rings first followed by the piston rings, and make sure that the ring join is not directly opposite the notch in the edge of the pump casing. See ⑨.

Mount the cylinders on the foot pedal by pressing the cross-pins over the hooks on the pedal and pushing towards the centre until they snap into place on the springs. See ⑩.

Place the cylinders on the pump chassis one at a time, pressing the first slightly down onto the piston before fitting the next. Make sure that the piston rings are correctly engaged all the way round. See ⑪.

Press the manifold pipe into position in the inlet opening, and then attach the bottom cover securely to the pump body. Finally, mount the suction tube and suction tip. See ⑫.

## 9. Test for correct function

- 1) Check that all four valve diaphragms and the manifold pipe in the container are mounted.
- 2) Block the suction tip with one finger. Depress one side of the pedal and check for vacuum by removing the finger.
- 3) Now block the suction tip once more and repeat process depressing the other side of the pedal. Check that a vacuum is built up.
- 4) Block the suction tip with one finger and pump several times. Wait approx. 10 seconds. When the finger is removed from the suction tip, the airflow being drawn into the pump should be audible.

## 10. Maintenance

All parts of the Ambu TwinPump are manufactured from high-strength, long-life materials. They therefore require no special maintenance or storage conditions.

# 11. Trouble-shooting

Defect	Possible cause	Correction
Suction diminished or non-existent and pedal is very easy to move, possibly on one side only.	Piston rings and/or O-rings not mounted or defective.	Mount or replace piston rings and/or O-rings.
	Valve diaphragms absent or defective.	Mount or replace valve diaphragms.

# 12. Spare parts - see ⑬

Cat. No.	Description	See
239 000 501	Ambu TwinPump: Pump chassis complete with valve diaphragms, O-rings and piston rings	13.1
240 000 501	Ambu TwinPump 1000: Pump chassis complete with valve diaphragms, O-rings and piston rings	13.1
000 239 010	Manifold pipe	13.2
000 239 005	Bottom cover	13.3
239 000 504	Cylinder with draw link and valve diaphragm	13.4
000 239 013	Valve diaphragm, separate	13.5
000 239 007	Piston ring	13.6
000 239 012	O-ring	13.7
239 000 508	Pedal with retaining springs	13.8
239 000 509	Suction tube, 135 cm, with combination suction tip and angle connector	13.9
239 000 510	Combination suction tip	13.10
000 239 006	Angle connector	13.11

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